



Medical and Laboratory Supplies

Since its founding in 1949, Heinz Herenz Medizinalbedarf GmbH has proven itself as innovative and renowned partner for customers in the health care sector in more than 90 countries all over the world. Market-orientated development of new innovations as well as constant advancement of our products, protected by own patents and trademark rights, enable us to supply goods of highest standards at attractive pricing.

We provide our customers with high quality CE-certified consumables and disposables for laboratories and hospitals, produced in accordance with the latest state of the art, and under consideration of the high ISO standards.

Extensive customer-oriented consultation and advice through our skilled team of experts is equally part of the conception of ourselves as well as punctual, prompt delivery just in time based on adequate stock-keeping.

Constant good quality of our products, in combination with fair market prices, have proved a recipe for the satisfaction of our customers. And this is not just a matter of good organization and our certified Quality Assurance System. Even more it is a matter of attitude.

We are deeply convinced to be able to meet your expectations and are looking forward to welcome you as our new customer.



Container

- Transport container for pathology
- Pathological vials
- Specimen container
- Sputum container, PS
- Faeces container
- Transport container
- Specimen container
- Universal cups
- Urine cups
- Snap-on cap

Heinz Herenz Medizinalbedarf GmbH

Rudorffweg 10
21031 Hamburg, Germany
Fon +49-40-739 20 40
Fax +49-40-730 41 48
Email: hhh@herenz.de
Web: www.herenz.de



Heinz Herenz Medizinalbedarf GmbH



Container

Transport container for pathology ▪ Pathological vials ▪ Specimen container ▪ Sputum container, PS
 Faeces container ▪ Transport container ▪ Specimen container ▪ Universal cups ▪ Urine cups ▪ Snap-on cap

Transport container for pathology



PP made, clear with screw cap (PE, not mounted), 35 ml container can be used as transport container for 20 ml, 100 ml container can be used as transport container for 50 ml. Available in: 20, 35, 50 and 100 ml

1131996-9

Pathological vials



Clear glass, washed, with aluminium screw cap, incl. black rubber seal. 7 ml, Bijou
 14 ml, McCartney
 28 ml, Universal Wide Neck

1070850-2

Specimen container



PP made, with screw cap, conical, autoclavable. approx. 40 ml, PP

1131819

Sputum container, PS



With pressure stopper PE, without spoon. 22 x 62 mm, 20 ml

1131000

Faeces container



PS With pressure stopper and plastic spoon at the lid., 22 x 63 mm. approx. 18 ml.

CE

1131001



PS 25 x 90 mm, with plastic spoon at the lid and screw cap. approx. 25 ml

CE

1131002



Universal container With mounted red screw cap. 30 ml, PP, with and without plastic spoon

CE

1131008

1131009



PP With spoon and screw cap PE, amber colour. 76 x 20 mm, 15 ml
 100 x 16 mm, 13 ml
 107 x 25 mm, 30 ml

CE

1131005-7

Transport container



For bacteriological samples, 110 x 30 x 30 mm, PE, with separate white screw cap, including absorbent paper.

Other colours of the cap by request.

1131859

Specimen container



Photo without stickers

Various types, PP made, each 60 ml, with screw cap (red), appr. 45 x 60 mm, with or without spoon and/or label

CE

1131003/4
 1131020-3

Universal cups



White, PS, environmentally friendly, 100 % recyclable. 150 and 200 ml

1131107/11

Urine cups



PP made, 125 ml, graduated to 100 ml, frosted field for writing, with screw cap (green/yellow), PE, sterile and nonsterile.

1131201-3

CE



PP made, 125 ml, graduated, unbreakable, frosted field for writing, nonsterile, for snap-on cap 1131121

CE

1131108

Snap-on cap



LDPE, red, with pour-out, easy to close, for urine cup 1131108.

CE

1131121