



## Medical and Laboratory Supplies

Since its founding in 1949, Heinz Herenz Medizinalbedarf GmbH has proven itself as innovative and renowned partner for customers in the health care sector in more than 90 countries all over the world. Market-orientated development of new innovations as well as constant advancement of our products, protected by own patents and trademark rights, enable us to supply goods of highest standards at attractive pricing.

We provide our customers with high quality CE-certified consumables and disposables for laboratories and hospitals, produced in accordance with the latest state of the art, and under consideration of the high ISO standards.

Extensive customer-oriented consultation and advice through our skilled team of experts is equally part of the conception of ourselves as well as punctual, prompt delivery just in time based on adequate stock-keeping.

Constant good quality of our products, in combination with fair market prices, have proved a recipe for the satisfaction of our customers. And this is not just a matter of good organization and our certified Quality Assurance System. Even more it is a matter of attitude.

We are deeply convinced to be able to meet your expectations and are looking forward to welcome you as our new customer.



### Heinz Herenz Medizinalbedarf GmbH

Rudorffweg 10  
21031 Hamburg, Germany  
Fon +49-40-739 20 40  
Fax +49-40-730 41 48  
Email: [hkh@herenz.de](mailto:hkh@herenz.de)  
Web: [www.herenz.de](http://www.herenz.de)



## Gloves

- Latex
- Vinyl
- Nitril
- Surgical
- PE
- Copolymer
- Wall mount rack

### Heinz Herenz Medizinalbedarf GmbH



# Gloves • Latex • Vinyl • Nitril • Surgical • PE • Copolymer • Wall mount rack

## Latex examination gloves



**powdered**  
Best quality, AQL 1.5.  
Size: XS / S / M / L / XL

1130483/4+1130499-1130501



**powderfree**  
Best quality, AQL 1.5.  
Size: S / M / L / XL

1130481/1130493-5

## Vinyl examination gloves



**powdered**  
Size: S / M / L / XL

1130480/1130497/1130498/113050



**powderfree**  
White  
Size: S / M / L / XL

1130539/1130490-2



**SUPERSTRETCH**  
**Powderfree**, more elastic than standard Vinyl gloves, AQL 1.5, EN 455.  
Size: XS / S / M / L / XL

1130570-4

## Nitril gloves



**blue**  
**powderfree**  
Size: XS / S / M / L / XL / XXL

1130530-4/1130550-2



**extra long execution**  
30 cm,  
Size: S / M / L

**nature**  
**powderfree**  
Size: S / M / L

1130535-7



**High Risk**  
Powderfree, 30 cm long, three times thicker than standard nitril gloves, ideal for handling of laboratory chemicals, exceptional protection against chemical exposure, AQL 0.65, EN 420, EN 374, EN 388, CE CAT 3.  
Size: S / M / L / XL / XXL / XXXL


1130575-80



Particularly suitable for technical trade, powderfree, AQL 1.5, EN 455. Colours: white, blue  
Size: S / M / L / XL

1130581-8

## Wall mount rack



Stainless steel, for glove dispenser boxes with customary dimensions.  
(max. 24 x 13 x 8,5 cm).

1130470



Translucent, acrylic glass,  
inner dimensions:  
250 x 135 x 95 mm.

1130471

## Gloves



**Neoprene gloves**  
**CE CAT III**, economic Neoprene-Latex composition, black, approx. 30 cm long, thick flocklining, for chemical and petrochemical applications, use in construction chemistry, paint shops et al.  
Size: 8-10

1110102/8



**Household gloves, latex**  
Suitable for foodstuff, flocklined, approx. 30 cm long.  
Colours: yellow, blue  
Size: S / M / L / XL

1130540-7



**Industrial gloves**  
Nitril, **CE-CAT III (3.0.0.1)**, green, 32 cms long, TOP-quality, anatomically shaped, structured inner surface, suitable for foodstuff handling, for use in the chemical industry and as protective glove.  
Size: XS / S / M / L / XL / XXL

1130589-94



**PE examination gloves**  
Embossed, polyethylene with 5 fingers, in 2 different packages.  
for men and women

1130496/1130503/1130518+19



**Surgical gloves**  
Latex made, sterile, powdered, packed in pairs.  
Size: 6,5 / 7 / 7,5 / 8 / 8,5

1130511-5



**Copolymer examination gloves**  
ETO-sterilized, non powdered, presented on bearing paper single packed.  
Size: S / M / L

1130505/1130509+10