



Medical and Laboratory Supplies

Since its founding in 1949, Heinz Herenz Medizinalbedarf GmbH has proven itself as innovative and renowned partner for customers in the health care sector in more than 90 countries all over the world. Market-orientated development of new innovations as well as constant advancement of our products, protected by own patents and trademark rights, enable us to supply goods of highest standards at attractive pricing.

We provide our customers with high quality CE-certified consumables and disposables for laboratories and hospitals, produced in accordance with the latest state of the art, and under consideration of the high ISO standards.

Extensive customer-oriented consultation and advice through our skilled team of experts is equally part of the conception of ourselves as well as punctual, prompt delivery just in time based on adequate stock-keeping.

Constant good quality of our products, in combination with fair market prices, have proved a recipe for the satisfaction of our customers. And this is not just a matter of good organization and our certified Quality Assurance System. Even more it is a matter of attitude.

We are deeply convinced to be able to meet your expectations and are looking forward to welcome you as our new customer.



Heinz Herenz Medizinalbedarf GmbH

Rudorffweg 10
21031 Hamburg, Germany
Fon +49-40-739 20 40
Fax +49-40-730 41 48
Email: hvh@herenz.de
Web: www.herenz.de



Gynaecology

- Gynobrush®
- Gynobrush® Set
- Ayres spatula
- Rovers® original
- Vaginal speculum
- Cyto-Fixx
- Sonosafe® probe covers
- Wooden applicator

Heinz Herenz Medizinalbedarf GmbH



Gynaecology

Gynobrush® • Gynobrush® Set • Ayres spatula • Rovers® original • Vaginal speculum • Cyto-Fixx
Ultrasound gel • Sonosafe® probe covers • Wooden applicator

Gynobrush®

suitable for gynaecological sample taking, non sterile



Gynobrush®

1032928



Gynobrush® plus
with soft plastic tip

1032929



Gynobrush® special
length approx. 20 cm, with specially
designed brush head

1032926



Gynobrush® combi
length approx. 24 cm combined
with a plastic made ayres spatula

1032927

Gynobrush® Set

Pap Smear Kit for gynaecology, non sterile, blister-packed



CE

1092929

each consisting of: 1 ayres spatula finger model F, 1 microscopic
slide (frosted field), 1 Gynobrush®, 1 cotton swab.

Ayres spatula



CE

1090703

Plastic (PP), length approx. 22 cm,
non sterile.



CE

1090701

**Made of wood, for gynaecology,
non sterile**
finger model F



CE

1090702

spatula model S

Brush suitable for gynaecological sample taking,

Rovers® original sterile and nonsterile

Rovers® Cervex-Brush®



1032924/5

CE

**Rovers® Cervex-Brush®
Combi**



1032940/1

CE

Rovers® Spatula



1032946/7

CE

Rovers® Endo Cervex-Brush®

with 600 thin hairs, sterile and non sterile



1032922/3

CE

**Rovers® Endo Cervex-
Brush® - S**

taking from narrow stenotic cervical os



1032942/3

CE

Rovers® Viba-Brush



1032944/5

CE

Vaginal speculum



CE

1131501-3
Disposable, sterilized by ethylene
oxide, colour coded.



1110294/5/7

acc. to Grave
Stainless steel 18/8,
glossy.

Cyto-Fixx Fixationsspray



Pumpspray for
fixation of cell
samples.

CE

1032921

Ultrasound gel transparent



Vacuum-manufactured gel
free from microbubbles and
acoustically precise for a wide
frequency range. Will not
corrode the ultrasonic trans-
ducer probes. Will not sensi-
tize or cause irritation to the
skin. Aldehyde-free.
No perfume, colour or UVsta-
bilizers added. For external
use only. Tested cytotoxicity
and dermatologically.

CE

1132050/6

Other packing sizes available

Sonosafe® probe covers

Covers for sonography sondes, without reservoir, latex made



CE

rolled, slightly powered,
Ø 34 mm, bulk packed

1130520



CE

siliconised or dry and slightly
powdered, rolled, single packed,
in aluminium foil, Ø 34 mm

1130521/3

Wooden applicator



WA CA 11

CE

1032318
1032418

With fixed cellulose bud DBGM,
made for CA-test, cellulose bud
approx. 35 x 5–6 mm.
Length: 25 cm
30 cm



WA CA 11 dünn

CE

1032519
1032518

With fixed cellulose bud DBGM,
made for CA-test, cellulose bud
approx. 35 x 4–5 mm.
Length: 25 cm
30 cm



WA 5

CE

1031618
1031718
1031818

Single cotton bud
(Ø approx. 5 mm, length
approx. 35 mm).
Length: 200 x 3 mm
250 x 3 mm
300 x 3 mm



WA 6

CE

1031918

Single cotton bud
(Ø approx. 4 - 5 mm, length
approx. 12 mm).
Length: 230 x 2,2 mm



WA 9

CE

1033001

Single cotton bud
(Ø approx. 4 - 5 mm, length
approx. 12 mm).
Length: 200 x 2,2 mm