



Medical and Laboratory Supplies

Since its founding in 1949, Heinz Herenz Medizinalbedarf GmbH has proven itself as innovative and renowned partner for customers in the health care sector in more than 90 countries all over the world. Market-orientated development of new innovations as well as constant advancement of our products, protected by own patents and trademark rights, enable us to supply goods of highest standards at attractive pricing.

We provide our customers with high quality CE-certified consumables and disposables for laboratories and hospitals, produced in accordance with the latest state of the art, and under consideration of the high ISO standards.

Extensive customer-oriented consultation and advice through our skilled team of experts is equally part of the conception of ourselves as well as punctual, prompt delivery just in time based on adequate stock-keeping.

Constant good quality of our products, in combination with fair market prices, have proved a recipe for the satisfaction of our customers. And this is not just a matter of good organization and our certified Quality Assurance System. Even more it is a matter of attitude.

We are deeply convinced to be able to meet your expectations and are looking forward to welcome you as our new customer.



OUR ASSORTMENT

 INSTRUMENTS	 STERISTOPPERS®	 COTTON-SWABS	 CARE	 CLOTHES
 LABORATORY	 REDUCTIONS	 TEST	 DISINFECTION	 GLOVES
 SPECIMEN CONTAINERS	 GLASS BOXES STAINING DISHES	 GLASS	 GYNAECOLOGY	 PIPETTE TIPS
 DISPOSABLES	 HAEMACYTOMETER	 PORCELAIN	 UROLOGY	 PIPETTES

Masks and glasses

- Face masks
- Filtering half masks
- Safety and leisure glasses
- Mask goggle
- Eye Shield
- Face Shield

Heinz Herenz Medizinalbedarf GmbH

Rudorffweg 10
21031 Hamburg, Germany
Fon +49-40-739 20 40
Fax +49-40-730 41 48
Email: h hh@herenz.de
Web: www.herenz.de



Heinz Herenz Medizinalbedarf GmbH



Masks and glasses

Face masks ▪ Filtering half masks ▪ Safety and leisure glasses
Mask goggle ▪ Eye Shield ▪ Face Shield

Face masks

CE



non woven material, breathable fabric, 3 ply, with 4 ties or rubber tie

Paper made, 1 ply, with rubber tie

1131974/5
1131816

Filtering half masks



without valve CE EN 149:2001, protection against solid and fluid particles. (FFP1 und FFP2)

CE

1131940/7



with valve CE EN 149:2001, protection against solid and fluid particles (FFP1 and FFP2) and protection against toxic particles (FFP3 only).

CE

1131941/2/5



flat-folded with valve (FFP2 and FFP3) and without valve (FFP1) CE EN 149:2001

CE

1131943/4/6

Safety and leisure glasses Arty

Polycarbonate lens with UV-filter 400, optical quality class 1.



1134100

Arty 250/F, CE EN 166 F, clear, fogfree



1134101

Arty 250/yellow, CE EN 1836, yellow, fogfree



1134102

Arty 250-small/F, CE EN 166, smaller model than Arty 250/F clear, fogfree



1134103

Arty 260/F, CE EN 166 F, clear, fogfree and scratchresistant



1134104

Arty 260/yellow, CE EN 166 F, yellow, fogfree and scratchresistant

Safety and leisure glasses Puma

Polycarbonate lens with UV-filter 400, optical quality class 1.



1134105

Puma 1262/orange, CE EN 166 F, orange



1134106

Puma 1262/smoke, CE EN 166 F, smoke



1134107

Puma 1262/yellow, CE EN 166 F, yellow

Safety and leisure glasses Sporty

Modern design, polycarbonate lens with UV-filter 400, interchangeable, optical quality class 1, fogfree, adjustable temple, soft and flexible temple ends.



1134108

Sporty/smoke CE EN 166F, smoke



1134109

Sporty/yellow CE EN 1836, yellow

Mask goggle



1134125

PANOLUX, CE EN 166 S, F, with indirect ventilation, flexible type with adjustable elastic band and special large field of vision, fogless and clear polycarbonate lens.

Eye Shield



1134150/1

Extremely lightweight and comfortable, latexfree, "wrap-around" design provides both front and side eye protection, distortion-free, opticalgrade lens, fits comfortably with a surgical mask and over prescription eyewear. EyeShield and frames are packed separately.

Face Shield



1134152/3

Extremely lightweight and comfortable, provides overhead-/front- and side-eye-protection, latexfree, adjustable, fog-resistant, fits comfortably with a surgical mask and over prescription eye-wear, includes a clear, snap-in nosepiece for support. FaceShield and frames are packed separately.