



## Medical and Laboratory Supplies

Since its founding in 1949, Heinz Herenz Medizinalbedarf GmbH has proven itself as innovative and renowned partner for customers in the health care sector in more than 90 countries all over the world. Market-orientated development of new innovations as well as constant advancement of our products, protected by own patents and trademark rights, enable us to supply goods of highest standards at attractive pricing.

We provide our customers with high quality CE-certified consumables and disposables for laboratories and hospitals, produced in accordance with the latest state of the art, and under consideration of the high ISO standards.

Extensive customer-oriented consultation and advice through our skilled team of experts is equally part of the conception of ourselves as well as punctual, prompt delivery just in time based on adequate stock-keeping.

Constant good quality of our products, in combination with fair market prices, have proved a recipe for the satisfaction of our customers. And this is not just a matter of good organization and our certified Quality Assurance System. Even more it is a matter of attitude.

We are deeply convinced to be able to meet your expectations and are looking forward to welcome you as our new customer.



### Heinz Herenz Medizinalbedarf GmbH

Rudorffweg 10  
21031 Hamburg, Germany  
Fon +49-40-739 20 40  
Fax +49-40-730 41 48  
Email: [hkh@herenz.de](mailto:hkh@herenz.de)  
Web: [www.herenz.de](http://www.herenz.de)



## Microscopic slides and Cover glasses

- cut / ground edges slides
- Superfrost slides
- economy slides
- diagnostic cell slides
- adhesive slides
- slides with cavity
- dispenser for slides
- borosilicate glasses
- economy types
- for haemocytometers
- circular

### Heinz Herenz Medizinalbedarf GmbH



# Microscopic slides and Cover glasses

cut / ground edges slides ▪ Superfrost ▪ economy slides ▪ diagnostic cell slides  
adhesive slides ▪ slides with cavity ▪ dispenser for slides ▪ borosilicate glasses  
economy types ▪ for haemacytometers ▪ circular

## Microscopic slides

cut or ground edges with or without frosted field



Washed and cleaned, 76 x 26 x 1 mm, in abrasion-proof boxes of 50 pcs.

1041000-1042100

CE

## Microscopic slides, Superfrost



Without cellophane, cut edges or ground edges. Colours: white, blue, pink, yellow, green, orange

1041001/3

CE

## Microscopic slides, economy

cut or polished edges 90°, with or without frosted field




cellophanized, approx. 76 x 26 x 1 mm.

1041050/1041150  
1042050/1042150

CE

## Diagnostic cell slides




1041014

100% wettable black colour, teflon (PTFE) coated slides 76 x 26 mm, with or without ground edges (depending on design), extremely hy-drophobic to prevent cross-contamination, resistant to acetone, methanol, PBS and most other chemical products, various standard designs. (special design and other colours available upon request, min. 2,000 pcs)

## Adhesive microscopic slides


### Superfrost Plus Gold



1043002

75 x 25 x 1 mm, extra white glass, ground 90°, white field for writing, packed 50 pcs, glass and box in cello.

### Superfrost Plus




1043001

75 x 25 x 1 mm, white field for writing, extra white glass, ground 90°.

Superfrost field for writing available upon request, colours: blue, pink, yellow, or green

### Polysine



1043000

75 x 25 x 1 mm, white field for writing, extra white glass, ground 90°.

## Dispenser for slides



1081603

ABS made, for max. 72 slides.

## Microscopic slides with cavity



1042001-3

Ground edges, 76 x 26 mm.  
1, 2 and 3 cavities

## Cover glasses non-ground



CE

Rectangular, No. 1 (0.13–0.16 mm thickness), **borosilicate glass**.



CE

Different sizes

## Cover glasses Economy non-ground



CE

No. 1 (0.13-0.16 mm thickness), **borosilicate glass**, also suitable for automated processes (cover slipper).



CE

No. 1 (0.13-0.16 mm thickness), **AR-glass**, for manual use

Different sizes

## Cover glasses optically plane



CE

1052604

Ground edges, thickness 0.4 mm, for haemacytometer.

## Cover glasses circular



1051199-1051211

Cut edges, No. 1.

Different sizes